



Item	External specification
IP type	Fuji Imaging Plate ST-VI 14×17"、14×14"、10×12"、8×10"、15X30cm
IP Cassette Type	IP Cassettetype CC/type LC
Connectable Modalities	FCRView V5.0 Medical and Veterinary
Network Connection	100BaseTX/10BaseT
Output to the Console	Pixel density 10pixels/mm
Processing Capacity	14 x 17" Approx. 17 IP/hr 14 x 14" Approx. 20 IP/hr 10 x 12" Approx. 21 IP/hr 8 x 10" Approx. 29 IP/hr 15 x 30cm Approx. 21 IP/hr NOTE : The time required for changing the cassette is assumed to be 0(zero) seconds. The time required for image erasure is assumed to be exposure performed at 25mR of X-ray dose.
Time required for IP feeding/loading	14 x 17" Approx. 211sec. 14 x 14" Approx. 180sec. 10 x 12" Approx. 172sec. 8 x 10" Approx. 125sec. 15 x 30cm Approx. 172sec. NOTE : The time required for image erasure is assumed to be exposure performed at 25mR of X-ray dose.
Time Required to Display Image	14 x 17" Approx. 52 sec. 14 x 14" Approx. 46 sec. 10 x 12" Approx. 44 sec. 8 x 10" Approx. 36 sec. 15 x30cm Approx. 44 sec.
External Dimensions WxDxH	23.6in × 15.7in × 30.7in
Weight	154lbs.
Power Consumption	170VA or less
Environment Condition (Operating conditions)	Temperature/Humidity ; 59°F(40-70%RH) to 86°F(15-70%RH) (no dew condensation) Atmospheric pressure ; 750-1060hPa
Environment Condition (Not operating)	Temperature: 0 - 113°F Humidity: 10 - 90%RH (no dew condensation) Atmospheric pressure: 500 - 1060hPa